

NBR1 (7H2) MOUSE MAB

Cat.#: N261353

Product Name: Anti-NBR1 (7H2) Mouse Monoclonal Antibody

Synonyms: 1A13B; MIG19; Nbr1

UNIPROT ID: Q14596

Background: Acts probably as a receptor for selective autophagosomal degradation of ubiquitinated targets.

Immunogen: Purified recombinant protein expressed in E.coli.

Applications: IHC-P

Recommended Dilutions: IHC: 1/50-1/100

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 7H2-4B2-6D6

MW: -

Isotype: IgG1

Purification: Affinity Purified

Species Reactivity: Human,Rat,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human Skeletal Muscle Tissue using NBR1 (7H2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
