

NDUFS6 RABBIT MAB

Cat.#: N262586

Product Name: Anti-NDUFS6 Rabbit Monoclonal Antibody

Synonyms: CI-13kA; CII3KDA; CI-13kD-A

UNIPROT ID: O75380

Background: Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.

Immunogen: Recombinant protein of human NDUFS6

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-9E1

MW: Calculated MW: 14 kDa; Observed MW: 14 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat,Hamster


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Tags & Cell Markers

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NDUFS6 in COS7, Raji, 3T3, Jurkat lysates using NDUFS6 antibody.



Western blot analysis of NDUFS6 in rat heart lysates using NDUFS6 antibody.