

NEUROMEDIN B RECEPTOR RABBIT PAB

Cat.#: N225532

Product Name: Anti-Neuromedin B Receptor Rabbit pAb

Synonyms: Neuromedin B receptor antibody; Bombesin Receptor 1 antibody; BBI antibody

UNIPROT ID: P28336

Background: Receptor for neuromedin-B.

Immunogen: Synthetic Peptide of NMBR

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 43 kDa; Observed MW: 50 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience yyy yyy Neuromedin B Receptor

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Neuromedin B Receptor in Human Brain, rat Brain lysates using NMBR(Bombesin Receptor) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
