

NFIB RABBIT MAB

Cat.#: N262603

Product Name: Anti-NFIB Rabbit Monoclonal Antibody

Synonyms: CTF; HMGIC/NFIB; NF-1/B; NF1-B; NFIB; NFIB2; NFI-B; NFI-RED; Nfib; NFIB2; NFIB3

UNIPROT ID: O00712

Background: Recognizes and binds the palindromic sequence 5'-TTGGC'NGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.

Immunogen: A synthetic peptide of human NFIB

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R03-0E5

MW: Calculated MW: 47 kDa; Observed MW: 47,55 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NFIB in C6, HeLa lysates using NFIB antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
