

## NGF RABBIT PAB

**Cat.#:** N225096

**Product Name:** Anti-NGF Rabbit pAb

**Synonyms:** NGF; NGFB; Beta-nerve growth factor; Beta-NGF

**UNIPROT ID:** P01138

**Background:** Nerve growth factor (NGF) is a small, secreted protein and member of the neurotrophin family of growth factors that promote neuronal cell survival and differentiation. Producing cells release NGF that bind and activate TrkA high affinity receptors to mediate NGF-driven signaling. NGF also binds to a low affinity p75 (NTR) receptors, which belong to the death receptor family.

**Immunogen:** The antiserum was produced against synthesized peptide derived from human NGF. AA range:33-82

**Applications:** WB,IHC-P,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 27 kDa; Observed MW: 27 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NGF in Jurkat lysates using NGF antibody.



Western blot analysis of NGF in various lysates using NGF antibody.