

NOXIN RABBIT PAB

Cat.#: N225194

Product Name: Anti-NOXIN Rabbit pAb

Synonyms: NOXIN; C11orf82; Nitric oxide-inducible gene protein

UNIPROT ID: Q8IXT1

Background: May be an anti-apoptotic protein involved in DNA repair or cell survival.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human NOXIN. AA range:111-160

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 112 kDa; Observed MW: 85 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NOXIN in 823, 293T, AD293, HeLa, HepG2 lysates using NOXIN antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
