

NR2C2AP RABBIT PAB

Cat.#: S216845

Product Name: Anti-NR2C2AP Rabbit Polyclonal Antibody

Synonyms: TRA16

UNIPROT ID: Q86WQ0 (Gene Accession - BC057837)

Background: May act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes.

Immunogen: Full length fusion protein

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human esophagus cancer tissue using 216845(NR2C2AP Antibody) at a dilution of 1/25(Nucleus).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
