

NTPCR RABBIT MAB

Cat.#: N262631

Product Name: Anti-NTPCR Rabbit Monoclonal Antibody

Synonyms: Clorf57; HCR-NTPase

UNIPROT ID: Q9BSD7

Background: Has nucleotide phosphatase activity towards ATP, GTP, CTP, TTP and UTP. Hydrolyzes nucleoside diphosphates with lower efficiency.

Immunogen: Recombinant protein of human Nucleoside triphosphate phosphohydrolase

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-6I7

MW: Calculated MW: 21 kDa; Observed MW: 21 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nucleoside triphosphate Phospho-hydrolase in K562, C6, 3T3, HeLa lysates using NTPCR antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
