

NUDT19 RABBIT MAB

Cat.#: N262632

Product Name: Anti-NUDT19 Rabbit Monoclonal Antibody

Synonyms: nudix (nucleoside diphosphate linked moiety X)-type motif 19; RP2

UNIPROT ID: A8MXV4

Background: Coenzyme A diphosphatase that mediates the hydrolysis of a wide range of CoA esters, including choloyl-CoA and branched-chain fatty-acyl-CoA esters. At low substrate concentrations medium and long-chain fatty-acyl-CoA esters are the primary substrates .

Immunogen: Recombinant protein of human NUDT19

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-9A5

MW: Calculated MW: 42 kDa; Observed MW: 42 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using NUDT19 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of NUDT19 (green) in HL-60 using NUDT19 antibody, and DAPI (blue).



Western blot analysis of NUDT19 in K562, rat Brain lysates using NUDT19 antibody.