

NUP107 RABBIT MAB

Cat.#: N263463

Product Name: Anti-NUP107 Rabbit Monoclonal Antibody

Synonyms: ODG6; NUP84; NPHS11; ODG6; GAMOS7

UNIPROT ID: P57740

Background: Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance (PubMed:12552102, PubMed:15229283). Required for the assembly of peripheral proteins into the NPC (PubMed:15229283, PubMed:12552102). May anchor NUP62 to the NPC (PubMed:15229283).

Immunogen: Recombinant protein of human Nup107

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-5F6

MW: Calculated MW: 106 kDa; Observed MW: 106 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Tags & Cell Markers

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nup107 in A549, HL-60, U251 lysates using NUP107 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
