

NUP62CL RABBIT MAB

Cat.#: N263138

Product Name: Anti-NUP62CL Rabbit Monoclonal Antibody

Synonyms: NUP62L; Nucleoporin-62 C-terminal-like protein

UNIPROT ID: Q9H1M0

Background: This gene encodes a protein containing a domain found in nucleoporins which are glycoproteins found in nuclear pore complexes. Alternative splicing results in multiple transcript variants.

Immunogen: A synthetic peptide of human NUP62CL

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-1A5

MW: Calculated MW: 21 kDa; Observed MW: 21 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NUP62CL in HeLa, Ramos lysates using NUP62CL antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
