

NUP88 RABBIT PAB

Cat.#: N225162

Product Name: Anti-NUP88 Rabbit pAb

Synonyms: NUP88; Nuclear pore complex protein Nup88; 88 kDa nucleoporin; Nucleoporin Nup88

UNIPROT ID: Q99567

Background: Component of nuclear pore complex.

Immunogen: Synthesized peptide derived from Nup88 . at AA range: 330-410

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 84 kDa; Observed MW: 84 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NUP88 in K562 lysates using Nup88 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
