

PARATHYROID HORMONE RABBIT PAB

Cat.#: N225837

Product Name: Anti-Parathyroid Hormone Rabbit pAb

Synonyms: hPTH; Parathormone; Parathyrin; Parathyroid hormone I; Prepro PTH; Preproparathyroid hormone; PTH1; PTH1 receptor; PTH1R; PTHR1

UNIPROT ID: P01270

Background: PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Stimulates [1-¹⁴C]-2-deoxy-D-glucose (2DG) transport and glycogen synthesis in osteoblastic cells.

Immunogen: A synthesized peptide derived from human Parathyroid Hormone

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 13 kDa; Observed MW: 36 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience Parathyroid Hormone

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Parathyroid Hormone in Human Parathyroid Hormone recombinant protein using Parathyroid Hormone antibody.