

PCGF6 RABBIT MAB

Cat.#: N262679

Product Name: Anti-PCGF6 Rabbit Monoclonal Antibody

Synonyms: RING finger protein 134; Mel18 and Bmi1-like RING finger

UNIPROT ID: Q99NA9

Background: Transcriptional repressor. May modulate the levels of histone H3K4Me3 by activating KDM5D histone demethylase. Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development.

Immunogen: Recombinant protein of mouse PCGF6

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-7J6

MW: Calculated MW: 40 kDa; Observed MW: 40 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Mouse, Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PCGF6 in PC-12 lysates using PCGF6 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
