

PDIA2 RABBIT PAB

Cat.#: N225774

Product Name: Anti-PDIA2 Rabbit pAb

Synonyms: P4HB; PDIA2; PDIR

UNIPROT ID: Q13087

Background: Acts as an intracellular estrogen-binding protein. May be involved in modulating cellular levels and biological functions of estrogens in the pancreas. May act as a chaperone that inhibits aggregation of misfolded proteins.

Immunogen: A synthesized peptide derived from human PDI

Applications: WB,FC

Recommended Dilutions: WB: 1/500-1/1000 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 58 kDa; Observed MW: 58 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PDI in HeLa lysates using PDIA2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
