

PFKFB3 RABBIT MAB

Cat.#: N262702

Product Name: Anti-PFKFB3 Rabbit Monoclonal Antibody

Synonyms: 6-bisphosphatase; 6-P2ase 3; Fructose-2; IPFK2; PFK/FBPase 3; PFK2; PFKFB3; uPFK 2

UNIPROT ID: Q16875

Background: Synthesis and degradation of fructose 2,6-bisphosphate.

Immunogen: Recombinant protein of human PFKFB3

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-7A4

MW: Calculated MW: 60 kDa; Observed MW: 60 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PFKFB3 in K562 lysates using PFKFB3 antibody.



Western blot analysis of PFKFB3 in Raw264.7, C6 lysates using PFKFB3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
