

PHLDA1 RABBIT MAB

Cat.#: N262710

Product Name: Anti-PHLDA1 Rabbit Monoclonal Antibody

Synonyms: PQ-rich protein; PQR protein; Proline- and histidine-rich protein; T-cell death-associated gene 51 protein

UNIPROT ID: Q8WV24

Background: Seems to be involved in regulation of apoptosis.

Immunogen: A synthetic peptide of human PHLDA1

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-9F4

MW: Calculated MW: 45 kDa; Observed MW: 45 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PHLDA1 in A375 lysates using PHLDA1 antibody.