

PHOSPHO-C MYB (SER11) RABBIT MAB

Cat.#: N263076

Product Name: Anti-Phospho-c Myb (Ser11) Rabbit Monoclonal Antibody

Synonyms: MYB; Transcriptional activator Myb; Proto-oncogene c-Myb

UNIPROT ID: P10242

Background: Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

Immunogen: A synthetic phosphopeptide corresponding to residues surrounding Ser11 of human v-Myb/c-Myb

Applications: WB,IHC-P,IP,ChIP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20 CHIP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-2D3

MW: Calculated MW: 72 kDa; Observed MW: 80 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Phosphorylated


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using vMyb/c-Myb (Phospho-Ser11) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Phospho-vMyb + c-Myb (Ser11) in Jurkat lysates using Phospho-c Myb (Ser11) antibody.