

PHOSPHO-ERK1/2 (THR202/TYR204)/(THR185/TYR187) RABBIT MAB

Cat.#: N262190

Product Name: Anti-Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187)
Rabbit Monoclonal Antibody

Synonyms: MAPK1/MAPK3

UNIPROT ID: P27361/P28482

Background: Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. MAPK1/ERK2 and MAPK3/ERK1 are the 2 MAPKs which play an important role in the MAPK/ERK cascade. They participate also in a signaling cascade initiated by activated KIT and KITLG/SCF. Depending on the cellular context, the MAPK/ERK cascade mediates diverse biological functions such as cell growth, adhesion, survival and differentiation through the regulation of transcription, translation, cytoskeletal rearrangements.

Immunogen: A synthetic phosphopeptide corresponding to residues surrounding Thr202/Tyr204 of human Erk1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-9E9

MW: Calculated MW: 44,42 kDa; Observed MW: 44,42 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Phosphorylated

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience, Prion disease

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Erk1 (pT202/pY204)/Erk2 (pT185/pY187) in C6, HeLa lysates using Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187) antibody.