

## PHOSPHO-RB2 P130 (SER952) RABBIT MAB

**Cat.#:** N262824

**Product Name:** Anti-Phospho-Rb2 p130 (Ser952) Rabbit Monoclonal Antibody

**Synonyms:** Rb2; P130

**UNIPROT ID:** Q08999

**Background:** Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation. Recruits and targets histone methyltransferases KMT5B and KMT5C, leading to epigenetic transcriptional repression. Controls histone H4 'Lys-20' trimethylation. Probably acts as a transcription repressor by recruiting chromatin-modifying enzymes to promoters. Potent inhibitor of E2F-mediated trans-activation, associates preferentially with E2F5. Binds to cyclins A and E. Binds to and may be involved in the transforming capacity of the adenovirus E1A protein. May act as a tumor suppressor.

**Immunogen:** A synthetic phosphopeptide corresponding to residues surrounding Ser952 of human Rb2 p130

**Applications:** WB,IHC-F,IHC-P,ICC/IF

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R02-5D6

**MW:** Calculated MW: 128 kDa; Observed MW: 128 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Phosphorylated


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Rb2 p130 (Phospho-Ser952) (green) in HeLa using Rb2 p130 (Phospho-Ser952) antibody, and DAPI (blue)



Western blot analysis of Rb2 p130 (Phospho-Ser952) in K562, HeLa lysates using Phospho-Rb2 p130 (Ser952) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Rb2 p130 (Phospho-Ser952) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.