

PHOSPHOLIPASE A2 ACTIVATOR PROTEIN RABBIT MAB

Cat.#: N262713

Product Name: Anti-Phospholipase A2 Activator Protein Rabbit Monoclonal Antibody

Synonyms: Glycerophosphatase; PLA2P; Plaa; PLAP

UNIPROT ID: Q9Y263

Background: Involved in the maintenance of ubiquitin levels.

Immunogen: A synthetic peptide of human Anti-Phospholipase A2 activator protein (PLAP)

Applications: WB,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-1D7

MW: Calculated MW: 87 kDa; Observed MW: 87 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospholipase A2 activator protein (PLAP) in 3T3, HeLa lysates using Phospholipase A2 activator protein (PLAP) antibody.



Immunocytochemistry analysis of Phospholipase A2 Activator Protein (green) in HeLa using Phospholipase A2 Activator Protein antibody, and DAPI (blue).