

PIAS1/2 RABBIT MAB

Cat.#: N263618

Product Name: Anti-PIAS1/2 Rabbit Monoclonal Antibody

Synonyms: PIAS1; PIAS2; PIAS3; Protein inhibitor of activated STAT protein 1

UNIPROT ID: O75925/O75928/Q9Y6X2

Background: O75928,Q9Y6X2. These proteins are involved in the pathway protein sumoylation, which is part of Protein modification.

Immunogen: A synthesized peptide derived from human PIAS1

Applications: WB, IHC-P, IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-7C1

MW: Calculated MW: 72 kDa; Observed MW: 76 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PIAS1+PIAS2 in HeLa, A549, HL-60, U2OS lysates using PIAS1+PIAS2 antibody.



Western blot analysis of PIAS1/PIAS2 in mouse testis, rat testis lysates using PIAS1/2 antibody.