

PKA C GAMMA RABBIT MAB

Cat.#: N262758

Product Name: Anti-PKA C gamma Rabbit Monoclonal Antibody

Synonyms: PRKACG; cAMP-dependent protein kinase catalytic subunit gamma; PKA C-gamma

UNIPROT ID: P22612

Background: Phosphorylates a large number of substrates in the cytoplasm and the nucleus.

Immunogen: Recombinant protein of human PRKACG

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-811

MW: Calculated MW: 40 kDa; Observed MW: 40 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PRKACG in MCF-7, C6, HeLa lysates using PRKACG antibody.



Western blot analysis of PRKACG in mouse testis, mouse lung lysates using PKA C gamma antibody.