

PLACENTAL LACTOGEN RABBIT PAB

Cat.#: N225211

Product Name: Anti-Placental Lactogen Rabbit pAb

Synonyms: chorionic somatomammotropin hormone 1; PL; CSA; CS-1; CSMT; hCS-1; hCS-A

UNIPROT ID: P0DML2

Background: Produced only during pregnancy and is involved in stimulating lactation, fetal growth and metabolism.

Immunogen: Synthetic peptide from human protein at AA range: 161-210

Applications: WB,IHC-P,ELISA

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 25 kDa; Observed MW: 25 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Placental Lactogen in mouse brain lysates using Placental Lactogen antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
