

PREB RABBIT MAB

Cat.#: N261621

Product Name: Anti-PREB Rabbit Monoclonal Antibody

Synonyms: Prolactin regulatory element-binding protein; Mammalian guanine nucleotide exchange factor mSec12

UNIPROT ID: Q9HCU5

Background: Guanine nucleotide exchange factor that specifically activates the small GTPase SAR1B.

Immunogen: Recombinant protein of human PREB

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R09-8B7

MW: Calculated MW: 45 kDa; Observed MW: 45 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PREB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PREB in K562, C6, 3T3, HeLa lysates using PREB antibody.