

PRMT3 RABBIT MAB

Cat.#: N262760

Product Name: Anti-PRMT3 Rabbit Monoclonal Antibody

Synonyms: HRMTIL3

UNIPROT ID: O60678

Background: Methylates (mono and asymmetric dimethylation) the guanidino nitrogens of arginyl residues in some proteins.

Immunogen: A synthetic peptide of human PRMT3

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-7K2

MW: Calculated MW: 60 kDa; Observed MW: 60 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PRMT3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PRMT3 in rat Brain, C6, 3T3, Hela lysates using PRMT3 antibody.