

RAB29 RABBIT MAB

Cat.#: N263519

Product Name: Anti-Rab29 Rabbit Monoclonal Antibody

Synonyms: RAB7L; RAB7L1

UNIPROT ID: O14966

Background: The small GTPases Rab are key regulators in vesicle trafficking.

Immunogen: Recombinant protein of human RAB29

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-3C0

MW: Calculated MW: 23 kDa; Observed MW: 23 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RAB29 in HeLa, A549, HL-60, U251 lysates using Rab29 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
