

RAIDD RABBIT MAB

Cat.#: N263523

Product Name: Anti-RAIDD Rabbit Monoclonal Antibody

Synonyms: CRADD; MGC9163; RAIDD; Death adaptor molecule RAIDD; Death domain containing protein CRADD

UNIPROT ID: P78560

Background: The receptor interacting protein RIP is a death domain-containing serine/threonine kinase which associates with FAS or the TNF-R1 binding protein TRADD. RAIDD (RIP-associated ICH-1/Ced-3 homologous protein with a death domain) has been identified as a RIP binding protein that also associates with members of the caspase family, providing a link between activation of the TNF-Rs and the triggering of the cysteine protease cascade. The amino-terminal domain of RAIDD shares significant homology with the prodomain of ICH-1 and mediates the binding of RAIDD to this cysteine protease.

Immunogen: A synthetic peptide of human RAIDD

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-1B3

MW: Calculated MW: 23 kDa; Observed MW: 23 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of RAIDD (green) in HeLa using RAIDD antibody, and DAPI (blue)



Western blot analysis of RAIDD in HeLa, A549, MCF-7, 293T lysates using RAIDD antibody.