

RANK RABBIT PAB

Cat.#: N225136

Product Name: Anti-RANK Rabbit pAb

Synonyms: CD265; FEO; LOH18CR1; mRANK; ODFR; OFE; OPTB7; OST5; PDB2; RANK; TNFRSF11A; TRANCER

UNIPROT ID: Q9Y6Q6

Background: Receptor for TNFSF11/RANKL/TRANCE/OPGL; essential for RANKL-mediated osteoclastogenesis. Involved in the regulation of interactions between T-cells and dendritic cells.

Immunogen: Synthetic peptide from human protein at AA range: 60-120

Applications: WB,IHC-P,ELISA

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 66 kDa; Observed MW: 66 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver using RANK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of RANK in Sodium Hela lysates using RANK antibody.