

SCN9A RABBIT PAB

Cat.#: S214653

Product Name: Anti-SCN9A Rabbit Polyclonal Antibody

Synonyms: PN1; ETHA; NENA; SFNP; FEB3B; NE-NA; GEFSP7; Nav1.7

UNIPROT ID: Q15858 (Gene Accession - NP_002968)

Background: This gene encodes a voltage-gated sodium channel which plays a significant role in nociception signaling. Mutations in this gene have been associated with primary erythralgia, channelopathy-associated insensitivity to pain, and paroxysmal extreme pain disorder.

Immunogen: Synthetic peptide of human SCN9A

Applications: ELISA, IHC

Recommended Dilutions: IHC: Oct-50; ELISA: 1000-2000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Metabolism, Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 214653(SCN9A Antibody) at a dilution of 1/20(Cell membrane).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the synthetic peptide and then with 214653(Anti-SCN9A Antibody) at dilution 1/20.