

TCTP RABBIT MAB

Cat.#: N262981

Product Name: Anti-TCTP Rabbit Monoclonal Antibody

Synonyms: Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23

UNIPROT ID: P13693

Background: Involved in calcium binding and microtubule stabilization.

Immunogen: A synthetic peptide of human TCTP

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-5D7

MW: Calculated MW: 20 kDa; Observed MW: 20 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of TCTP (green) in HeLa using TCTP antibody, and DAPI (blue)



Western blot analysis of TCTP in K562, rat Brain, C6, 3T3, HeLa lysates using TCTP antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
