

TENSIN 2 RABBIT PAB

Cat.#: N225197

Product Name: Anti-Tensin 2 Rabbit pAb

Synonyms: TENC1; KIAA1075; TNS2; Tensin-like C1 domain-containing phosphatase; C1 domain-containing phosphatase and tensin homolog; C1-TEN; Tensin-2

UNIPROT ID: Q63HR2

Background: Tyrosine-protein phosphatase which regulates cell motility proliferation and muscle-response to insulin.

Immunogen: The antiserum was produced against synthesized peptide derived from human Tensin-2. AA range:439-488

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 153 kDa; Observed MW: 125 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tensin 2 in A549 3T3 lysates using Tensin 2 antibody.