

TFG RABBIT MAB

Cat.#: N263027

Product Name: Anti-TFG Rabbit Monoclonal Antibody

Synonyms: TF6; HMSNP; SPG57; TRKT3

UNIPROT ID: Q92734

Background: Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

Immunogen: A synthetic peptide of human TRK fused gene

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-8F5

MW: Calculated MW: 43 kDa; Observed MW: 43 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TRK fused gene antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TRK fused gene in C6, 3T3 lysates using TFG antibody.



Immunocytochemistry analysis of TRK fused gene (green) in HeLa using TRK fused gene antibody, and DAPI (blue)