

TFG RABBIT PAB

Cat.#: N225732

Product Name: Anti-TFG Rabbit pAb

Synonyms: TF6; HMSNP; SPG57; TRKT3

UNIPROT ID: Q92734

Background: Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

Immunogen: A synthesized peptide derived from human TRK fused gene

Applications: WB,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 43 kDa; Observed MW: 55 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TRK fused gene in HeLa lysates using TFG antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
