

TGF ALPHA RABBIT PAB

Cat.#: N225944

Product Name: Anti-TGF alpha Rabbit pAb

Synonyms: ETGF; TFGA; TGF A; TGF-alpha; Tgfa; Wa1; Waved 1

UNIPROT ID: P01135

Background: TGF alpha is a mitogenic polypeptide that is able to bind to the EGF receptor and to act synergistically with TGF beta to promote anchorage-independent cell proliferation in soft agar.

Immunogen: A synthesized peptide derived from human TGF alpha

Applications: WB,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 17 kDa; Observed MW: 17 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TGF alpha in (1) HepG2 lysates; (2) Human fetal liver lysates using TGF alpha antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
