

THAP12 RABBIT PAB

Cat.#: S218810

Product Name: Anti-THAP12 Rabbit Polyclonal Antibody

Synonyms: DAP4; THAP0; PRKRIR; P52rIPK

UNIPROT ID: O43422 (Gene Accession - BC021992)

Background: Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.

Immunogen: Fusion protein of human THAP12

Applications: ELISA, IHC

Recommended Dilutions: IHC: 40-200; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human ovarian cancer tissue using 218810(THAP12 Antibody) at a dilution of 1/45(Nucleus).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
