

THIOREDOXIN DOMAIN CONTAINING PROTEIN 17 RABBIT PAB

Cat.#: N225171

Product Name: Anti-Thioredoxin Domain Containing Protein 17 Rabbit pAb

Synonyms: TXNDC17; TXNL5; Thioredoxin domain-containing protein 17; 14 kDa thioredoxin-related protein; TRP14; Protein 42-9-9; Thioredoxin-like protein 5

UNIPROT ID: Q9BRA2

Background: Disulfide reductase. May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions. Modulates TNF- α signaling and NF- κ -B activation. Has peroxidase activity and may contribute to the elimination of cellular hydrogen peroxide.

Immunogen: Synthesized peptide derived from TRP14 . at AA range: 20-100

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 14 kDa; Observed MW: 14 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Thioredoxin Domain Containing Protein 17 in HL-60, Raji, K562, SW480, A549, 22RV1 lysates using Thioredoxin Domain Containing Protein 17 antibody.