

THO COMPLEX SUBUNIT 1 RABBIT PAB

Cat.#: N226049

Product Name: Anti-THO Complex Subunit 1 Rabbit pAb

Synonyms: Death domain containing protein p84N5; HPR1; hTREX84; Nuclear matrix protein p84; p84N5; THO complex 1; Tho1; Thoc1

UNIPROT ID: Q96FV9

Background: Regulates transcriptional elongation of a subset of genes. Participates in an apoptotic pathway which is characterized by activation of caspase-6, increases in the expression of BAK1 and BCL2L1 and activation of NF-kappa-B.

Immunogen: A synthetic peptide of human Nuclear Matrix Protein p84

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 76 kDa; Observed MW: 84 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Tags & Cell Markers

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nuclear Matrix Protein p84 in 293T, PC-12, 3T3 lysates using THO Complex Subunit 1 antibody.