

TMA16 RABBIT PAB

Cat.#: S218832

Product Name: Anti-TMA16 Rabbit Polyclonal Antibody

Synonyms: C4orf43

UNIPROT ID: Q96EY4 (Gene Accession - BC011842)

Background: The function of this protein remains unknown.

Immunogen: Fusion protein of human TMA16

Applications: ELISA, IHC

Recommended Dilutions: IHC: 100-200; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human esophagus cancer tissue using 218832(TMA16 Antibody) at a dilution of 1/105(Nucleus and Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human esophagus cancer tissue is first treated with the fusion protein and then with 218832(Anti-TMA16 Antibody) at dilution 1/105.