

TNFA(GOLIMUMAB) MAB

目录: 28580

产品名称: TNF α (Golimumab) Monoclonal Antibody

基因符号: DIF,TNF- α ,TNFA,TNFSF2,cachexin,cachectin,TNF α

描述: TNF α (golimumab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: ELISA, Flow Cyt

推荐稀释度: ELISA 1:5000-10000, Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human TNF α

成分: Lyophilized from sterile PBS, pH 7.4. 5% - 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 μ g/mL TNF α (golimumab biosimilar) mAb (28580) on Expi293 cells transfected with Human TNF α protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).




Figure 2. ELISA plate pre-coated by 2 μ g/mL (100 μ l/well) Human TNF α ± Protein, His Tag(11199) can bind TNF α (golimumab biosimilar) mAb(28580) in a linear range of 3.20-80 ng/mL.