

TNFRSF10C RABBIT PAB

Cat.#: S216468

Product Name: Anti-TNFRSF10C Rabbit Polyclonal Antibody

Synonyms: LIT; DCRI; TRID; CD263; TRAILR3; TRAIL-R3; DCRI-TNFR

UNIPROT ID: O14798 (Gene Accession - BC125041)

Background: The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular TRAIL-binding domain and a transmembrane domain, but no cytoplasmic death domain. This receptor is not capable of inducing apoptosis, and is thought to function as an antagonistic receptor that protects cells from TRAIL-induced apoptosis. This gene was found to be a p53-regulated DNA damage-inducible gene. The expression of this gene was detected in many normal tissues but not in most cancer cell lines, which may explain the specific sensitivity of cancer cells to the apoptosis-inducing activity of TRAIL.

Immunogen: Fusion protein of human TNFRSF10C

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cancer, Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver cancer tissue using 216468(TNFRSF10C Antibody) at a dilution of 1/30(Cytoplasm and Cell membrane).



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using 216468(Anti-TNFRSF10C Antibody) at a dilution of 1/30.