

TOLL-LIKE RECEPTOR 13 RABBIT PAB

Cat.#: N225530

Product Name: Anti-Toll-Like Receptor 13 Rabbit pAb

Synonyms: Toll-Like Receptor 13; Toll like receptor 13 antibody

UNIPROT ID: Q6R5N8

Background: Component of innate and adaptive immunity that recognizes and binds 23S rRNA from bacteria. TLRs (Toll-like receptors) control host immune response against pathogens through recognition of molecular patterns specific to microorganisms. Acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response.

Immunogen: Recombinant Protein of TLR13

Applications: IHC-P

Recommended Dilutions: IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: -

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded mouse Spleen Tissue using Toll-Like Receptor 13 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.