

TRANSMEMBRANE PROTEIN 161B RABBIT MAB

Cat.#: N262998

Product Name: Anti-Transmembrane Protein 161B Rabbit Monoclonal Antibody

Synonyms: Transmembrane protein 161B

UNIPROT ID: Q8NDZ6

Background: Essential for maintaining normal cardiac rhythm in the developing heart and for neonatal survival (By similarity).

Immunogen: Recombinant protein of human TMEM161B

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-6J4

MW: Calculated MW: 55 kDa; Observed MW: 55 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TMEM161B in K562, HeLa lysates using Transmembrane Protein 161B antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
