

TRANSMEMBRANE PROTEIN 72 RABBIT PAB

Cat.#: N225244

Product Name: Anti-Transmembrane Protein 72 Rabbit pAb

Synonyms: Transmembrane protein 72; Kidney-specific secretory protein of 37 kDa; C10orf127; KSP37

UNIPROT ID: A0PK05

Background: This gene encodes a transmembrane protein which may be expressed specifically in the kidney.

Immunogen: Synthesized peptide derived from human Tmem72
Polyclonal

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 30 kDa; Observed MW: 30 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Transmembrane Protein 72 in various lysates using Transmembrane Protein 72 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
