

TRERF1 RABBIT PAB

Cat.#: S214825

Product Name: Anti-TRERF1 Rabbit Polyclonal Antibody

Synonyms: RAPA; BCAR2; TREP132; TReP-132; HSA277276; dJ139D8.5

UNIPROT ID: Q96PN7 (Gene Accession - NP_277037)

Background: This gene encodes a zinc-finger transcriptional regulating protein which interacts with CBP/p300 to regulate the human gene CYP11A1. Alternative splicing results in multiple transcript variants encoding different isoforms.

Immunogen: Synthetic peptide of human TRERF1

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cancer

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human esophagus cancer tissue using 214825 (TRERF1 Antibody) at a dilution of 1/20 (Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human esophagus cancer tissue is first treated with the synthetic peptide and then with 214825 (Anti-TRERF1 Antibody) at dilution 1/20.