

TRIM15 RABBIT PAB

Cat.#: S219476

Product Name: Anti-TRIM15 Rabbit Polyclonal Antibody

Synonyms: RNF93; ZNFB7; ZNF178

UNIPROT ID: Q9C019 (Gene Accession - BC038585)

Background: The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to the cytoplasm. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.

Immunogen: Fusion protein of human TRIM15

Applications: ELISA, IHC

Recommended Dilutions: IHC: 100-300; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling, Cell Biology


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human esophagus cancer tissue using 219476 (TRIM15 Antibody) at a dilution of 1/90 (Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human esophagus cancer tissue is first treated with the fusion protein and then with 219476 (Anti-TRIM15 Antibody) at dilution 1/90.



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using 219476 (Anti-TRIM15 Antibody) at a dilution of 1/90.



In comparison with the IHC on the left, the same paraffin-embedded Human tonsil tissue is first treated with fusion protein and then with D227296 (Anti-TRIM15 Antibody) at dilution 1/90.