

## TRMT11 RABBIT PAB

**Cat.#:** N225202

**Product Name:** Anti-TRMT11 Rabbit pAb

**Synonyms:** TRMT11; C6orf75; MDS024; tRNA; guanine(10)-N2)-methyltransferase homolog; tRNA guanosine-2'-O-methyltransferase TRMT11 homolog

**UNIPROT ID:** Q7Z4G4

**Background:** Catalytic subunit of an S-adenosyl-L-methionine-dependent tRNA methyltransferase complex that mediates the methylation of the guanosine nucleotide at position 10 (m<sup>2</sup>G<sub>10</sub>) in tRNAs.

**Immunogen:** Synthesized peptide derived from the C-terminal region of human TRMT11.

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 53 kDa; Observed MW: 53 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** -

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TRMT11 in HepG2 , JK, Colo, Hela lysates using TRMT11 antibody.