

UBN2 RABBIT PAB

Cat.#: S221574

Product Name: Anti-UBN2 Rabbit Polyclonal Antibody

Synonyms:

UNIPROT ID: Q6ZU65 (Gene Accession - NP_775840)

Background: Ubinuclein 2 (UBN2) is a human ortholog of the yeast protein Hpc2p. It possesses an Hpc2 (histone promoter control 2)-related domain which is highly conserved.

Immunogen: Synthetic peptide of human UBN2

Applications: ELISA, IHC

Recommended Dilutions: IHC: 20-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human tonsil tissue using 221574(UBN2 Antibody) at a dilution of 1/20(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human tonsil tissue is first treated with the synthetic peptide and then with 221574(Anti-UBN2 Antibody) at dilution 1/20.