

UFC1 RABBIT MAB

Cat.#: N263052

Product Name: Anti-UFC1 Rabbit Monoclonal Antibody

Synonyms: NEDSG; HSPC155

UNIPROT ID: Q9Y3C8

Background: E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.

Immunogen: Recombinant protein of human UFC1

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-5E3

MW: Calculated MW: 19 kDa; Observed MW: 19 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of UFC1 (green) in HeLa using UFC1 antibody, and DAPI (blue)



Western blot analysis of UFC1 in K562, HeLa lysates using UFC1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
