

## UFC1 RABBIT MAB

**Cat.#:** N263051

**Product Name:** Anti-UFC1 Rabbit Monoclonal Antibody

**Synonyms:** NEDSG; HSPC155

**UNIPROT ID:** Q9Y3C8

**Background:** E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.

**Immunogen:** Recombinant protein of human UFC1

**Applications:** WB,IHC-P

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R08-5E1

**MW:** Calculated MW: 19 kDa; Observed MW: 19 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of UFC1 in K562, rat Brain, C6, 3T3 lysates using UFC1 antibody.